



# Flowpad-Flowstretch

## A generic toolkit for microfluidics and point-of-care diagnostics

### What is Flowpad-Flowstretch?

FlowPad features disposable, credit card format microfluidic cartridges with fluidic channels designed to perform a selection of tests outside of a traditional lab environment. The integrated protocol can be complex by including all sample preparation steps—concentration, lysis, purification—as well as bio-analysis and detection procedures—qPCR, RPA, LAMP and ELISA.

An instrument leverages a chip carrier, a functional implementation system based on the microfluidic cartridge design. This instrument integrates a valve, a pump and a mixing chamber and is embedded into a toolkit. This toolkit delivers the diagnosis. It integrates actuators, a Peltier heater, magnetic actuators and on-board reagent recovery, as well as an optical sensing system, and associated software programs.

Both instrument and microfluidic cartridges are custom-designed based on clinical or manufacturing requirements.

### Applications

FlowPad finds applications in human medicine, agri-agro and environment industries:

- Sample preparation at the sampling point  
—peptide, siRNA isolation
- Point-of-care testing and diagnostics  
—bacteria, viruses, biomarkers and proteins.
- Companion diagnostic
- Cell encapsulation—spheroids, cells, etc.
- Assessment of Chemical, Biological, Radiological and Nuclear (CBRN) risks.

## What's new?

CEA-Leti's Flowpad-Flowstretch technology opens a new path to the future of testing:

- Tests in less than 2 hours: from raw biological sample to results
- Versatile, easy-to-use toolbox
- Compatible with multiple complex protocols—biological, chemical, etc.—and materials—silicon, glass, COC, polymers, etc.
- Reduces development and prototyping time
- Easily scalable for high-volume production at low-cost. CEA-Leti has perfected cartridge design and small-scale production.
- Full protocol integration

Flowpad-Flowstretch's Microfluidic cartridge:

- Are user-friendly features plug-and-play fluidic connections with easy and efficient flow control
- Operates with valves and a pump integrated into the cartridge
- Are compatible with injection molding



## Publications

- **Quantitative biological assays with on-chip calibration using versatile architecture and collapsible chambers.** C. Parent, F. Boizot, M. Cubizolles, N. Verplanck, J.-L. Achard, Y. Fouillet
- **Microfluidic device integrating a network of hyper-elastic valves for automated glucose stimulation and insulin secretion collection from a single pancreatic islet.** C. Quintard, E. Tubbs, J.-L. Achard, F. Navarro, X. Gidrol, Y. Fouillet
- **Peps: An Innovative Microfluidic Device for Bedside Whole Blood Processing before Plasma Proteomics Analyses.** B. Gilquin, M. Cubizolles, R. Den Dulk, F. Revol-Cavalier, M. Alessio, C.-E. Goujon, C. Echampard, G. Arrizabalaga, A. Adrait, M. Louwagie, P. Laurent, F. Navarro, Y. Couté, M.-L. Cosnier, V. Brun



## What's next?

CEA-Leti is currently working in scaling to high-volume manufacturing both the cartridge and the instrument with industrial partners. New application: multiplex measurement card.

## Interested in this technology?

Contact:  
**Valérie Roux**  
[valerie.roux@cea.fr](mailto:valerie.roux@cea.fr)  
 +33 438 780 426

**CEA-Leti, technology research institute**

17 avenue des Martyrs, 38054 Grenoble Cedex 9, France

[cea-leti.com](http://cea-leti.com)

   @CEA-Leti

